

Design and performance of the X-123 compact X-ray and Gamma-ray spectroscopy system

R. Redus, A. Huber, J. Pantazis, T. Pantazis, D. Sperry

Amptek, Inc 14 DeAngelo Dr, Bedford MA 01730

This Web presentation is an extended version of the paper presented at the 2006 Room Temperature Semiconductor Detector Workshop, at the IEEE Nuclear Science Symposium



Outline

Design

- System
- Detector and preamp
- Digital processor
- Power supplies & packaging

Performance

- Resolution
- High Count Rates
- Stability
- Spectroscopy

Applications

X123 Design

Compact integrated system

- Detector and preamplifier
 Detector thermoelectrically cooled
- Digital pulse processor with shaping amplifier and MCA
- All power supplies
- Small size
 7 x 10 x 2.5 cm³ (2.7 x 3.9 x 1 in³)
- Low power (1.2W)
- Light weight (180 g or 6.3 oz)
- Simple to operate
 - USB communications
 - +5VDC







X123 Design



Block diagram of X123 (and related systems)

- Complete spectroscopy system: detector, preamp, digital processor (with shaping, MCA, and logic), comm ports, power supplies, packaging
- Same building blocks as Amptek's DP4 and PX4
 - Animation shows relationship between Amptek's digital processing products
- Detector options include Si-PIN, CdTe, and CdTe-stack

X123 Design

X123, PX4, DP4

- All are complete spectroscopy systems
- All share detectors, digital processing, other core technologies
- Targeted at different applications

PX4 DIGITAL PULSE PROCESSOF

XR-100CR

X-RAY DETECTOR



X123 for compact, packaged system

TEK

AMP

®



DP4 Digital Pulse Processor Pulse Shaping + MCA Power Supply PA-210 Preamplifier AXR, X-Ray Detector

PX4 and XR100 for benchtop & laboratory

Thermoelectrically Cooled Solid State Detector

- Reasons for thermoelectric cooling
 - Reduces shot noise and thermal noise
 - For CdTe, allows higher bias
 - Cooling invisible to user
- Cooler achieves >80°C differential
 - Now using 2 stage coolers in all products
 - 215K for lab use
 - 240K for field use (at ambient of 45°C)
- FET and feedback components on cooler
 - Leakage currents as low as 5 fA
 - Low stray capacitance, reduced EMI pickup



AMP

®

AMP TEK

X123 Detectors

Three detectors for different energy ranges

- Low energies (X-rays)
 - Si-PIN detector
 - Recommended from 1.5 to 25-40 keV
 - Resolution: 150 eV FWHM at 5.9 keV
- Medium energies
 - CdTe diode
 - Recommended from 25-40 to 100 keV
 - Resolution: 600 eV FWHM at 59.5 keV
- Gamma-Rays
 - CdTe diode stack
 - Recommended from 100 to 600 keV
 - Resolution: 6 keV at 662 keV
- All based on thermoelectric cooling
- All are compatible with X123



Plot compares a ⁵⁷Co spectrum from Si-PIN and CdTe detectors: the Si-PIN has much better resolution but above 30 keV much lower sensitivity.

Low energy X-rays

- Detector : Si-PIN
- Thickness: 200 to 1000 μ m
- Area: 5 to 25 mm²
- Feedback: Capacitive reset into the HV node
- Be window: 0.5 to 4 mil
- Internal multilayer collimator
 Prevents interactions near guard rings
- Best resolution 149 eV FWHM ENC = 100 eV FWHM \Rightarrow 12 electron rms Measured with 6 mm² x 500 µm detector at 215K and 20 µsec shaping time



AMP

®

Med energy X-rays

- Detector : CdTe Schottky diode
- Thickness: 750 to 1000 μ m
- Area: 9 to 25 mm²
- Feedback: Current divider
- Bias: 300 to 1000V
- Photofraction >90% to 80 keV
- Resolution
 - ENC as low as 500 eV FWHM
 - 750 ev FWHM at 59.5 keV
 - 1200 eV FWHM at 122 keV



AMP

®

Gamma-rays

- Detector : CdTe stack
 - Multiple planar elements
 - Summed into single preamp
 - Sensitivity arises from entire volume Thermoelectric Cooler
 - Max charge transport length 0.75 mm
 - Operates like single detector
- Thickness: 2.25 or 3.75 mm
- Area: 25 mm²
- Feedback: Resistive
- Bias: 1 to 1.4 kV

<1% FWHM at 662 keV

AMP

Page 10 of 40

®

Plot compares intrinsic efficiency of the detectors. At low energies, the Be window limits sensitivy. The Si-PIN have good efficiency to 30 keV, CdTe to 100 keV, and the CdTe-stack to 500 keV.

Recent Detector Developments

- Cooler achieves >80°C differential
 - 2 stage coolers in all products
 - 215K for lab use, 240K for field use (at ambient of 45°C)
- Thicker active depth for Si-PIN
 - 1000 $\mu m,$ fully depleted, with areas of 7 and 25 mm^2
 - Resolution 185 eV FWHM
- Si-PIN with internal multilayer collimator
 - Collimator prevents interactions near guard rings and outer edge
 - Multilayer design prevents characteristic X-rays contaminating signal

Work in Progress

- Si drift diodes
 - Will increase count rates to reduce measurement time
- Reset preamp for CdTe
 - Will combine lowest noise with good count rate stability

FPGA pipeline logic

- Integrated MCA & pulse selection logic
- Same core on DP4, PX4, and X123

Advantages of digital processing

- Flexible
 - Many shaping times: 24 log steps from 0.8 to 102 μ sec peaking
 - Adjustable settings for BLR, RTD, etc
 - All adjustments are in s/w, via serial connection
- Faster than analog
 - No dead time for pulse digitization
 - Finite Impulse Response improves baseline stability
 - Symmetric shape enhances throughput
- Pulse shape closer to optimum for signal to noise
- Improved stability and reproducibility

Trapezoidal Pulse Shape

- No long tail so reduced pile-up
- Finite impulse response reduces baseline shifts
- Symmetry reduces PUR interval
 Dead time 1.25 τ_{peak}+τ_{flat}
- Improved ENC for same pulse duration (FWHM)

Fast Channel

- 600 nsec resolving time
- Used for
 - Pile-up rejection,
 - Risetime discrimination
 - Measurement of incoming count rate

rval

Digital (Trapezoidal)

Analog (Quasi-Triangular)

Digital processing features

- Gain
 - Coarse gain range x11 to x110 (4 ranges)
 PX4 has larger coarse gain range, fully overlapping
 - Fine gain range 0.75 to 1.25, 10 bit resolution
 - All gain changes software selectable over serial bus
- System Conversion Gain
 - For Si-PIN at 1k channels, 7 to 110 eV/channel \Rightarrow 7 to 110 keV full scale
 - For CdTe, 20 to 310 eV/channel \Rightarrow 20 to 310 keV full scale
 - For CdTe stack, 50 to 3100 eV/channel \Rightarrow 50 to 3100 keV full scale
- Baseline Restoration
 - Asymmetric algorithm stabilizes negative noise peaks
 - Software selectable parameters (64 settings) to tune for application

MultiChannel Analyzer Features

- MCA Performance
 - Commandable to 256, 512, 1024, 2048, 4096k, 8192 channels.
 - 3 bytes (24 bits) per channel 16.7M counts
 - Minimum Acquisition Time < 10 msec
 - Data transfer time for 1k channels in 10 msec (USB) or 500 msec (RS-232)
 - Conversion time: None, part of digital pipeline

- Acquisition Modes

- Standard MCA Mode: Integrated pulse height spectrum
- Delta Mode: Pulse height spectrum from most recent 1 sec acquisition
- Repetitive Mode: After each acquisition, saves data to file and then restarts.
- MCS Mode: Each channel has number of counts in a time bin, >10 msec per bin

– Other features

- Presets: Time, total counts, counts in an ROI, counts in a channel
- External controls: Gate and buffer select inputs

Command and Control

- Software control over serial interfaces
 - Sends commands, starts and stops acquisition
 - Reads out and displays spectra and count rates
- Protocols supported
 - RS232 and USB standard
 - Can be tailored for Ethernet
 - I²C and μ P pins for direct h/w interfaces
- Software
 - Amptek's standard ADMCA software for data acquisition and control
 - Interface to with analysis packages (XRF-FP available from Amptek)
 - Library (DLL) of drivers for custom software

AMP TEK

X123 Outputs

- Primary output is the spectrum
 - Packet of spectral data
 - Typically 1k or 2k channels, updated 1 Hz
- Secondary output is oscilloscope display
 - Captured waveform sent over serial interface
 - Computer displays waveform, for setup and debugging
- Other outputs available on digital processor board
 - Standard outputs on DP4, PX4. X123 requires custom connector
 - Analog outputs: shaped pulse, ADC input, fast pulse, test pulser
 - Digital outputs: Timing trigger, reset flag, pile-up flag, peak flag, etc
 - Eight single channel analyzer outputs (LLD, ULD set in software)

X123 Power Design

Architecture

- Multiple switching supplies
- Low voltages: +3.3V digital, +/-5V analog, +/-8V preamp

Standard, commercial PWM supplies

High voltage: +100 to +1500VDC at microamps

Switch mode supply with Cockroft-Walton

- Cooler supply: Closed loop regulation of temperature

Commercial PWM supply with temperature feedback

Configurations

- X123: Output voltages set by internal jumpers and pots
- PX4: Voltages set by serial command

X123 Power Design

Input

- +5VDC nominal
- Actual range is 4V to 6V for battery operations
- Typical steady state 1.2W (240 mA at 5V)
- First minutes, 1.8W (360 mA) until detector is cooled

Power Budget

- 0.4W for Peltier cooler (typical steady state at 25°C)
- 0.4W for digital processor
- 0.2W for analog circuits
- 0.2W for power supplies

X123 Packaging

Packaging is critical!

- Charge sensitive preamp, with noise of 12 e- rms, located with cm of 0.3A switching supply
- Must remove 0.5W from detector with minimal $\Delta T_{ambient}$
- Machined AI box for EMI, thermal
- Careful grounding, layout, and shielding
- Connections
 - USB, serial, power standard
 - Others available
- Mounting kit shown

AMP

®

Si-PIN ⁵⁵Fe spectra

Best resolution 149 eV @ 5.9 keV

CdTe spectra

- Resolution 0.75 keV @ 59.5 keV

CdTe stack spectra

Resolution 5 keV @ 662 keV

Si-PIN Noise Components

- Typical results with a 13mm² detector
- Series noise dominates below 25 μ sec peaking time
- Parallel noise unimportant
- Fano broadening dominates >10 keV

CdTe Noise Components

- Typical results with a 25 mm² x0.75 mm detector
- Noise corner near 6.4 μsec peaking time
- Noise dominates below 30-50 keV
- Hole tailing dominates above this

What determines the energy resolution?

$$\Delta E^{2} = F \bullet E + ENC^{2} + F_{tail}(E) = \left(F \bullet E\right) + \left(\frac{n_{series}}{\tau} + n_{1/f} + n_{parallel}\tau\right) + F_{tail}(E)$$

- Add in quadrature: ΔE is square root of sum of squares, not linear sum.
- Statistical fluctuations
 - Contribution scales as \sqrt{FE} , F is Fano factor, about 0.1
 - A function of energy only. Dominates at high energy for Si-PIN
- Electronic noise
 - Independent of energy so dominates at low energy
 - A function of the shaping time, with a minimum at the "noise corner" where the series and parallel terms are equal.
 - Series noise is proportional to total input capacitance and dominates at short shaping times.
 - Parallel noise arises from detector leakage current and dominates at long shaping times.
 - 1/f (pink) noise does not vary with shaping time.
- Hole tailing is important for CdTe only, not Si-PIN
 - Complicated function of energy and voltage and varies from one detector to next
 - Dominates at high energies. Improved with increased bias voltage

X123 Sensitivity

Si-PIN Sensitivity

- For Si-PIN detectors, sensitivity vs. resolution trade-off
- 6 mm² x 500 μ m best resolution
- 25 mm² x 500 μm large area
- 13 mm² x 1000 μm has sensitivity to higher energies
- Must be matched to application

Multi-Element Fluorescence with 500 μm and 1000 μm thick Si-PIN Detector

Page 27 of 40

X123 High Count Rates

Dead Time and Throughput Different from Analog

- In analog, "dead time" specified is for peak digitization.
- In digital, there is no dead time for this, only for pulse overlap \Rightarrow much shorter dead time and so higher throughput
- In digital, pulse is symmetric and timing well defined \Rightarrow shorter and better known dead time intervals.
- In analog, fast signal used to detect pile-up but fast rate not recorded. In digital, directly measure the fast rate \Rightarrow more accurate correction.

Throughput

- Much better than analog system for same τ_{peak}
- Measured results matches equations very well
- Fast channel permits direct dead time correction

$$R_{Out} = R_{In} \left(e^{-R_{In} \left(\tau_{Peak} + \tau_{Flat} \right)} \right) \qquad \text{PUR Off}$$

$$R_{Out} = R_{In} \left(e^{-(2.375)(R_{In})(\tau_{Peak} + \tau_{Flar})} \right) \qquad \text{PUR On}$$

$$R_{In} = \frac{R_{FAST}}{1 - R_{Fast} \tau_{Fast}} \quad \text{Dead Time Correction}$$

X123 High Count Rates

Spectral performance at high rates

- We recommend operation <50% dead time
- Baseline shifts <0.2% to 50% dead time
- Resolution changes similarly
- Spectra here are up to 70% dead time
- Double peaks visible to 95% dead time

2006 Symposium on Room Temp. Detectors

AMP

Page 30 of 40

®

X123 Stability

Drift over time

- Top plot shows X123-Si data
- Bottom plot shows X123-CdTe
- Gain fluctuations consistent with temperature coefficient
 - 60 ppm/°C for Si-PIN
 - 30 ppm/°C for CdTe
- Count rate follows decay
- PX4 data taken for 1 month
- No long term drift observed
- Slope within std dev of zero

AMP

2006 Symposium on Room Temp. Detectors

®

X-Ray Fluorescence

- Method used to non-destructively identify and measure which elements are in a sample
- Physics:
 - Incoming X-ray or γ-ray ejects inner shell electron from an atom
 - Excited atom decays to ground state
 - During decay, atom emits characteristic X-ray with energy equal to difference in shells
 - Each element has a unique pattern of transitions so a unique, characteristic set of X-rays
- Presence of X-rays ⇒ presence of element
 Need sufficient energy resolution to separate peaks
- Rate of X-rays ⇒ quantity of sample
 Need accurate measurement of counts
 Need to correct for attenuation, scattering, etc

X-Ray Fluorescence: Sample Spectra

- Typical result for lead
- 88 keV incident photons excite sample
 - Some scattered photons reach the detector
 - Lead K-shall has 88 keV binding energy, so photoelectrons are ejected
 - Decays from L shells (13-15 keV) lead to characteristic X-ray emissions at 72, 75, and 85 keV
 - Vacancies in L shells are filled by transitions from M shells
 - Generates characteristic X-rays at 10, 12, and 15 keV
- No damage to sample
- Penetrates beyond surface layers
- Rapid and accurate

Sample Spectrum

- Used Amptek XR100, not X123
- Taken by Mars Pathfinder
- Spectrum shows peaks present in sample and relative areas
- Yields elemental composition of Martian rocks

AMP

2006 Symposium on Room Temp. Detectors

®

ROHS Verification

(Reduction of Hazardous Substances)

- EU requires commercial electronics to be free of lead, cadmium, mercury, etc. as of 7/06
- OEMs need to verify compliance of vendors and suppliers
- X123 Si-PIN can reliably quantify the critical elements
- Compact packaging, easy connections provide simple, COTS solution for end users

AMP

2006 Symposium on Room Temp. Detectors

®

Art and Archeology

- XR100 currently used for XRF analysis
- X123 provides same performance in much more compact form

AMP TEK

XRF of Michelangelo's David

Page 36 of 40

Industrial process control

- XRF provides key quality measurements
- X123 permits simple, integrated solution for end users

AMP

2006 Symposium on Room Temp. Detectors

Page 37 of 40

®

X-Ray Tube Quality Assurance

- Stable tube output is very important for many applications
- Mammography particular demanding
- X123 CdTe provides high sensitivity over the energy range of most tubes
- Validate beam at start of day
- Track beam flux over time
- Track spectral stability over time
- Compact, integrated X123 makes integration into system easy

AMP

R

Radioistope identification

- X123 CdTe-Stack
- Good energy resolution to 750 keV
- Complements higher sensitivity detectors such as NaI(TI)
- Small size, low power, simple interface simplify system integration

AMP

2006 Symposium on Room Temp. Detectors

®

X123

Conclusions

- Successfully integrated Peltier cooled detectors, digital pulse processing, and power supplies into a compact package
- Works well for Si-PIN, and for 1 mm CdTe, and for CdTe stacks to cover a wide energy range
- Provides high resolution of analog benchtop systems
- Digital processor provides much higher throughput and stability than analog benchtop system
- Compact size, low power, and simple interface simplify system integration for the end user
- For more information, go to www.amptek.com